


BIOPSY NEEDLES

Ultracore Biopsy Needle - Manual

Ultracore is used for performing core biopsies of soft tissues.



- Ultrasonically visible tip facilitates precise placement.
- CM marking for proper positioning.
- Wide guage and length range.

QTY/Box: 5

Order No.	G	Length (cm)
UCB 14/3 UCB 14/4.5 UCB 14/6	14	3" 4.5" 6"
UCB 16/3 UCB 16/4.5 UCB 16/6	16	3" 4.5" 6"
UCB 18/3 UCB 18/4.5 UCB 18/6	18	3" 4.5" 6"



SURGIMEDIK