


SPECIAL RADIOLOGY

Small Bowel Enema Catheter Open / Closed Tip

Used for double contrast display of small intestine.



- Round tip for atraumatic introduction.
- Soft radiopaque, flexible catheter for quick introduction.
- Heavy duty torque wire with manipulation knob.

Small Bowel Enema Catheter Closed Tip

QTY/Box: 1

Order No.	Description
SBEC-C 12/6-G	SBEC, Closed Tip, 12F, 6 sides holes with Torque wire, 150 cm
SBEC-C 14/6-G	SBEC, Closed Tip, 14F, 6 sides holes with Torque wire, 150 cm
SBEC-C 12/6	SBEC, Closed Tip, 12F, 6 sides holes 150 cm
SBEC-C 14/6	SBEC, Closed Tip, 14F, 6 sides holes 150 cm

Small Bowel Enema Catheter Open Tip

QTY/Box: 1

Order No.	Description
SBEC-O 12/0-G	SBEC, Open Tip, 12F, 0 sides holes with Torque wire, 150 cm
SBEC-O 14/0-G	SBEC, Open Tip, 14F, 0 sides holes with Torque wire, 150 cm
SBEC-O 12/0	SBEC, Open Tip, 12F, 0 sides holes 150 cm
SBEC-O 14/0	SBEC, Open Tip, 14F, 0 sides holes 150 cm

