

SNARES

Polysnare

Used for the endoscopic excision of polypectomy in conjunction with a monopolar electro-surgical unit.



- Wider range of design for physician choice
- PTFE sheath
- Wire braided design for easy maneuvering
- 1:1 torque for handle wire movement
- Banana Plug Connector

QTY /Box: 1

Snare Type	5 Fr		7 Fr	
	160 cm	240 cm	160 cm	240 cm
Small Oval	PS-SO-5/160	PS-SO-5/240	PS-SO-7/160	PS-SO-7/240
Medium Oval	PS-MO-5/160	PS-MO-5/240	PS-MO-7/160	PS-MO-7/240
Large Oval	PS-LO-5/160	PS-LO-5/240	PS-LO-7/160	PS-LO-7/240
Crescent	PS-C-5/160	PS-C-5/240	PS-C-7/160	PS-C-7/240
Hexagonal	PS-HX-5/160	PS-HX-5/240	PS-HX-7/160	PS-HX-7/240