

STONE EXTRACTION BALLOON

SEB Stone Extraction Balloon

Used for removal of biliary stones and to facilitate injection of contrast medium while occluding the duct with the balloon.



- Double lumen feature provides access to the guidewire to the desired location during extraction and also facilitates guidewire exchange.
- Heavy duty, symmetrical and low-profile balloons.

QTY /Box: 1

| Order No. | Guidewire Compatibility | Balloon | | Length (cm) | Fr |
|-----------|-------------------------|---------------|---------------|-------------|----|
| | | Capacity (ml) | Diameter (mm) | | |
| SEB7/200 | 0.035 | 1 | 8.5 | 200 | 7 |
| | | 2 | 12 | | |
| | | 4 | 15 | | |
| | | 5 | 18 | | |

Caution:

This product contains natural rubber latex which may cause allergic reaction in some individuals.