


### R+ Ureteral Stents (Open Tip, Closed Tip, Intra operative)

Used to temporarily facilitate drainage of urine from kidney to bladder by trained physician.



- A Full new line of ureteral stent with enhanced radiopacity for easier visualization & patient comfort
- 4.7 Fr Stents with long pusher for use via Ureteroscope

#### Stent + Pusher

QTY. / Box : 10 / 25

Fr	Length (cm)	GW	Open Taper Tip	Closed Taper Tip	Intra operative (both ends closed taper)
4.7	16	0.035"	SOT 4.7/16 - R+	SS 4.7/16 - R+	SSI 4.7/16 - R+
	24		SOT 4.7/24 - R+	SS 4.7/24 - R+	SSI 4.7/24 - R+
	26		SOT 4.7/26 - R+	SS 4.7/26 - R+	SSI 4.7/26 - R+
	28		SOT 4.7/28 - R+	SS 4.7/28 - R+	SSI 4.7/28 - R+
	30		SOT 4.7/30 - R+	SS 4.7/30 - R+	SSI 4.7/30 - R+
6	16	0.038"	SOT 6/16 - R+	SS 6/16 - R+	SSI 6/16 - R+
	24		SOT 6/24 - R+	SS 6/24 - R+	SSI 6/24 - R+
	26		SOT 6/26 - R+	SS 6/26 - R+	SSI 6/26 - R+
	28		SOT 6/28 - R+	SS 6/28 - R+	SSI 6/28 - R+
	30		SOT 6/30 - R+	SS 6/30 - R+	SSI 6/30 - R+
7	16	0.038"	SOT 7/16 - R+	SS 7/16 - R+	SSI 7/16 - R+
	24		SOT 7/24 - R+	SS 7/24 - R+	SSI 7/24 - R+
	26		SOT 7/26 - R+	SS 7/26 - R+	SSI 7/26 - R+
	28		SOT 7/28 - R+	SS 7/28 - R+	SSI 7/28 - R+
	30		SOT 7/30 - R+	SS 7/30 - R+	SSI 7/30 - R+

#### Combinations :

With Guidewire : Add Suffix - G

With Suture : Add Suffix - S

With Clamp : Add Suffix - C

With Hydrophilic Coating : Add Suffix - H

#### e.g. :

Stent + Pusher + Guidewire : SS 6/26 - R+ - G

Stent + Pusher + Suture : SS 6/26 - R+ - S

Stent + Pusher + Clamp : SS 6/26 - R+ - C

Stent + Pusher + Suture + Guidewire + Clamp : SS 6/26 - R+ - S - G - C

#### Domestic:

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.

Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

#### International:

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com