

Endopyelotomy Stents

The Endopyelotomy Stent is used for temporary internal drainage from the UPJ to Bladder following incision of a stricture.



- Soft and flexible for patient comfort
- Biocompatible polysurethane with low encrustation tendency and superior radiopacity for enhanced visualization
- Large diameter at the incised UPJ provides for optimal healing
- 13 cm extended segment

QTY. / Box : 5

Taper Size (Fr)		Legth Between	GW	Order No.
From	To			
6	12	24	0.038"	EP 6/12-24 EP 6/12-26 EP 6/12-28
		26		
		28		
7	14	24	0.038"	EP 7/14-24 EP 7/14-26 EP 7/14-28
		26		
		28		

* **Combinations** : (With Guidewire : Add Suffix - G)

* **Each stent is supplied with pusher**

Domestic:

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.
Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

International:

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com