

Percutaneous Nephrostomy Drainage Set-Pigtail

Used for percutaneous placement of a pigtail catheter in the pelvis for nephrostomy drainage.



- A convenient complete set for percutaneous placement of drainage catheter-pigtail design
- Material provides superior biocompatibility and drainage
- Highly radiopaque for precise placement
- Locking loop design available for better anchoring
- Softer polyurethane
- Contact us for other sizes

Set includes :

- Initial Puncture Needle - 3 Part,
- J Tip fixed core PTFE Coated 0.038" Guidewire
- Fascial Dilators
- Pigtail Drainage Catheter 30 cm
- Urine Bag Connector
- Scalpel
- Catheter Fixing Disk
- Stop Cock

QTY. / Box : 1

Fr	Order No.
7	PDS-P7
8	PDS-P8
9	PDS-P9
10	PDS-P10
12	PDS-P12
14	PDS-P14

Combinations :

(With Locking Loop System : Add Suffix : LS)

(With Lunderquist Guidewire instead of PTFE Guidewire : Add Suffix - L)

(With Hydrophilic Coating : Add Suffix : H)

Domestic:

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.

Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

International:

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com