


### Ultima Percutaneous Nephrostomy Drainage Sets

Used for percutaneous placement of a pigtail catheter in the pelvis for nephrostomy drainage.



- A convenient complete set for percutaneous placement of drainage catheter-pigtail design
- Material provides superior biocompatibility and drainage
- Highly radiopaque for precise placement
- Locking loop design available for better anchoring
- Firm polyurethane
- Contact us for other sizes

Set includes :

- Ultima Drainage Catheter
- Initial Puncture Needle - 3 Part,
- J Tip fixed core PTEF Coated 0.038" Guidewire
- Fascial Dilators
- Urine Bag Connector
- Scalpel
- Catheter Fixing Disk
- Stop Cock

QTY. / Box : 1

Fr	Order No.
8	U-PDS-P8/30
9	U-PDS-P9/30
10	U-PDS-P10/30
12	U-PDS-P12/30
14	U-PDS-P14/30

**Combinations :**

(With Locking Loop System : Add Suffix : LS)

(With Lunderquist Guidewire instead of PTFE Guidewire : Add Suffix - L)

(With Hydrophilic Coating : Add Suffix : H)

**Domestic:**

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.

Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

**International:**

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com