


Percutaneous Nephrostomy Drainage Set-Malecot

Used for percutaneous placement of a malecot catheter in renal for nephrostomy drainage.



- A convenient complete set for percutaneous placement of drainage catheter-malecot design (ideally suitable for the drainage of thick, viscous fluids.
- Material provides superior biocompatibility and drainage
- Highly radiopaque for precise placement
- 30 cm length
- Soft polyurethane
- Contact us for other sizes

Set includes :

- Malecot Drainage Catheter 30 cm
- Initial Puncture Needle - 3 Part,
- J Tip fixed core PTFE Coated 0.038" Guidewire
- Fascial Dilators
- Urine Bag Connector
- Scalpel
- Catheter Fixing Disc
- Stop Cock

QTY. / Box : 1

| Fr | Order No. |
|----|-----------|
| 10 | PDS-M10 |
| 12 | PDS-M12 |
| 14 | PDS-M14 |
| 16 | PDS-M16 |

Combinations :

(With Underquist Guidewire instead of PTFE Guidewire : Add Suffix - L)

(With Hydrophilic Coating : Add Suffix : H)

Domestic:

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.

Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

International:

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com