

Suprapubic Puncture & Drainage Set

Used for percutaneous bladder drainage.



- Ideal when drainage is required in emergency situation
- Helpful for sexually active patients or when Urethral catheterisation is not an option
- Comprehensive procedural kit for ease of use
- Available in select places only

- Suprapubic Puncture Set includes : Suprapubic Trocar, Sleeve, Scalpel
- Suprapubic Drainage Set includes : Suprapubic Trocar, Sleeve, Scalpel, Forey Balloon Catheter


QTY. / Box : 10

QTY. / Box : 5

Fr	Order No. (Suprapubic Puncture Set)	Order No. (Suprapubic Drainage Set with Silicone Elastomer)
8	SPS 8	SDS 8-SE
12	SPS 12	SDS 12-SE
16	SPS 16	SDS 16-SE

Suprapubic Metal Tip Puncture & Drainage Set

Used for nephrostomy drainage through an existing track.



- Suprapubic Puncture Set includes :
Suprapubic Trocar Metal Tip, Sleeve, Scalpel
- Suprapubic Drainage Set includes :
Suprapubic Trocar Metal Tip, Sleeve, Scalpel, Foley Balloon Catheter

QTY. / Box : 10

QTY. / Box : 5

Fr	Order No. (Suprapubic Puncture Set)	Order No. (Suprapubic Drainage Set with Silicone Elastomer)
12	SPS-12-M	SDS-12-SE-M
16	SPS-16-M	SDS-16-SE-M

Domestic:

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.

Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

International:

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com