

UreTrack Ureteral Balloon Dilators

The Uretrack Ureteral Balloon Dilators are recommended for the dilatation of the urinary track.



- Accepts 0.035" Guidewire
- 5 Fr Shaft
- 18 atm rated burst pressure
- 5 mm (15 Fr) dilation
- Supplied with Stop Cock
- Option of set with inflation Device 20cc

QTY. / Box : 1

Order No.	Order No. (with Inflation Device)	Balloon Inflated OD	Balloon Length	Shaft Fr	Length (cm)
UBD 5/4 UBD 5/10	UBD 5/4-IDA UBD 5/10-IDA	5mm (15 Fr) 5mm (15 Fr)	4 10	5	70

TU Transureteroscopic Balloon Dilator

The Transureteroscopic Ureteral Balloon Dilators are recommended for the dilatation of the ureter through the 5 Fr-working channel of an Ureteroscope.



- Provided with stylet
- Accepts 0.018" Guidewire
- 18 atm rated burst
- 4 mm (12 Fr) X 4 cm balloon, 3.5 Fr Shaft diameter
- Supplied with Stop Cock
- Option of set with Inflation Device 20 cc

QTY. / Box : 1

Order No.	Order No. (With Inflation Device)	Balloon Inflated OD	Balloon Length	Shaft Fr	Length (cm)
TU 4/4	TU 4/4-IDA	4 mm (12 Fr)	4	3.5	70

Domestic:

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.

Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

International:

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com