

## Inflation Device Advanced

Used with ballon dilation catheters to monitor inflation pressure.




- Ergonomic design facilitates excellent eye to hand co-ordination
- Efficient intra cam locking mechanism allows to maintain pressure and guarantee fast and automatic release
- Superior design with grained surface and slip free lever ensures maximum grip even with wet hands
- Pressure gauge ensures fatigue free readability
- Upto 30 / 26 atm Analogue pressure range
- Registers immediately without pressure decay or drift


QTY. / Box : 1

Order no.	Description
IDA 20	Inflation Device Advanced 20 cc

## Urinary Tract Occlusion Balloon Catheter

Used for UPJ Occlusion to prevent fragment migration during percutaneous nephrostomy procedures.



- Facilitates hydrodilatation of the collecting system, to aid in percutaneous puncture and increase internal manipulating space
- Supplied with stylet and stop cock
-  Contain natural rubber Latex which may cause allergy in some individuals

QTY. / Box : 1

Fr	Length (cm)	Order No.
3	65	UTOB 3/65

**Domestic:**

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.

Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

**International:**

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com