


Nitiflex Nitinol Stone Basket - Flatwire

The Stone Basket has been designed for endoscopic manipulation, dislodging, engaging and removal of renal, ureteral stones.




- Made of Nitinol, Flatwire design - 4 wires
- Available in 2.4 F, 3 F, 4 F sizes to meet all requirements

QTY. / Box : 1

Fr	Length (cm)	Order No.	Order No. (with Polyimide Sheath)
2.4	90	NITI 4W/2.4F-90-F	NITI 4W/2.4F-90-F-PA
	120	NITI 4W/2.4F-120-F	NITI 4W/2.4F-120-F-PA
3	90	NITI 4W/3F-90-F	NITI 4W/3F-90-F-PA
	120	NITI 4W/3F-120-F	NITI 4W/3F-120-F-PA
4	90	NITI 4W/4F-90-F	NITI 4W/4F-90-F-PA
	120	NITI 4W/4F-120-F	NITI 4W/4F-120-F-PA

Nitiflex Nitinol Stone Basket - Helical

The Stone Basket has been designed for endoscopic manipulation, dislodging, engaging and removal of renal, ureteral stones.



- made of Nitinol, Helical design - 4 wires
- Available in 2.4 F, 3 F, 4 F sizes to meet all requirements

QTY. / Box : 1

Fr	Length (cm)	Order No.	Order No. (with Polyimide Sheath)
2.4	90	NITI 4W/2.4F-90-H	NITI 4W/2.4F-90-H-PA
	120	NITI 4W/2.4F-120-H	NITI 4W/2.4F-120-H-PA
3	90	NITI 4W/3F-90-H	NITI 4W/3F-90-H-PA
	120	NITI 4W/3F-120-H	NITI 4W/3F-120-H-PA
4	90	NITI 4W/4F-90-H	NITI 4W/4F-90-H-PA
	120	NITI 4W/4F-120-H	NITI 4W/4F-120-H-PA

Domestic:

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.

Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

International:

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com