


## Trapper No Tip Nitinol Stone Basket

The Stone Basket has been designed for endoscopic manipulation, dislodging, engaging removal of renal, ureteral stones.



- Nitinol construction of shaft & basket
- No tip design
- Dismantleable handle
- Available in selected places only
- 1.8 Fr Stone Basket is with Nitinol Shaft having flexible portion for use with flexible Ureteroscope

QTY. / Box : 1

| Fr  | Length (cm) | Order No.                         | Order No.<br>(with Polyimide Sheath)    |
|-----|-------------|-----------------------------------|---|
| 1.8 | 90<br>120   | -                                 | TNB 4W/1.8F-90-PA<br>TNB 4W/1.8F-120-PA |
| 2.4 | 90<br>120   | TNB 4W/2.4F-90<br>TNB 4W/2.4F-120 | TNB 4W/2.4F-90-PA<br>TNB 4W/2.4F-120-PA |
| 4   | 90<br>120   | TNB 4W/3F-90<br>TNB 4W/3F-120     | TNB 4W3F-90-PA<br>TNB 4W3F-120-PA       |
| 4   | 90<br>120   | TNB 4W/4F-90<br>TNB 4W/4F-120     | TNB 4W/4F-90-PA<br>TNB 4W/4F-120-PA     |

**Domestic:**

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.  
Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

**International:**

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com