


Supremo Nitinol Hydrophilic Guidewires

Used to access ureter during routine and difficult cases.




- Nitinol Non-Kink core
- Hydrophilic coated
- in straight and Angled design
- Super maneuverability

QTY. / Box : 1

Tip	Diameter Tip	Coated	Length (cm)	Order No.
Straight	0.025"	Hydrophilic	150	SHSG 150-25
	0.032"			SHSG 150-32
	0.035"			SHSG 150-35
Angled	0.032"			SHAG 150-32
	0.035"			SHAG 150-35

Cobra Nitinol Guidewire

Used for percutaneous entry into vessel using the seldinger technique or endoscopic access.



- Non - Kink Nitinol Shaft
- Black and White stripes for easier visualization
- Straight flexible tip

QTY. / Box : 1

Length (cm)	Diameter	Order No.
150	0.025"	CGS 150-25
	0.032"	CGS 150-32
	0.035"	CGS 150-35

Domestic:

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.

Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

International:

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com