


## Lunderquist Guidewires

For percutaneous entry into vessel using the seldinger technique or for endoscopic access.




- Firm and Easy to introduce
- Shapeable tip facilitates accurate placement even in difficult anatomy
- Available in 100 cm length
- 3 mm J-Tip

QTY. / Box : 1

Tip	Diameter Tip	Length (cm)	Order No.
3 mm J	0.035"	100	LQJ 100-35-3

## HyGuide Hydrophilic Guidewires

Used for percutaneous entry into vessel using the seldinger technique or for endoscopic access.



- Exacting tolerances of each wire facilitates easy reliable instrument delivery
- Hydrophilic coating and bonded polymer Jacket permit smooth navigation through tortuous vasculature
- Non-Kink Nitinol Shaft

QTY. / Box : 5

Description	Order No.
Straight Fixed Core, Hydrophilic Coated 150 cm - 0.035"	HGSF 150-35

### Domestic:

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.  
Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

### International:

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com