

Schuller Guidewires

For percutaneous entry into vessel using the seldinger technique or for endoscopic access.



- Firm and Easy to introduce
- Shapeable tip facilitates accurate placement even in difficult anatomy
- Available in 100 cm length
- 3 mm J-Tip

QTY. / Box : 1

Tip	Floppy Part	Rigid Shaft	Length (cm)	Shaft	Order No.
J-Curved , 3 mm	60 cm	30 cm	100	0.035"	SHJ 100-30-60

Standard Straight & J Tip Guidewires

Used for percutaneous entry into vessel using the seldinger technique or for endoscopic access.

Straight - TIP


QTY. / Box : 5



Description	Order No.
Straight Fixed Core, Teflon Coated, 150 cm - 0.018"	SFT 150-18
Straight Fixed Core, Teflon Coated, 150 cm - 0.025"	SFT 150-25
Straight Fixed Core, Teflon Coated, 150 cm - 0.028"	SFT 150-28
Straight Fixed Core, Teflon Coated, 150 cm - 0.032"	SFT 150-32
Straight Fixed Core, Teflon Coated, 150 cm - 0.035"	SFT 150-35
Straight Fixed Core, Teflon Coated, 150 cm - 0.038"	SFT 150-38

J - TIP

QTY. / Box : 5



Description	Order No.
J Tip, 3 mm, Fixed Core, Teflon Coated, 150 cm, 0.032"	JFT 150-32-3
J Tip, 3 mm, Fixed Core, Teflon Coated, 150 cm, 0.035"	JFT 150-52-3
J Tip, 3 mm, Fixed Core, Teflon Coated, 150 cm, 0.038"	JFT 150-38-3

Domestic:

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.
Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

International:

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com