

Flexi Cath

1


Straight Dual Lumen Catheter

- Soft tip for easier insertion
- Bodysoft polyurethane material
- One step “Stepper Dilator” included
- Available in catheter and kit version
- Seldinger technique insertion


Colour coded luers



Secured anchoring



Soft Tip



Printed Priming Volume



Bodysoft polyurethane material



Catheter Contents :

1 ea Flexi-Cath, 1 ea Stepper Dilator 12 F, 2 ea Injection Caps.

Kit Contents :

1 ea Flexi-Cath, 1 ea Stepper Dilator 12 F, 2 ea Injection Caps.

1 ea Introducer Needle 18/7, 1 ea Guidewire JFS 80-38, 1 ea 5ml Syringe, 1 ea Scalpel.

Order No.	Description	Qty. / Box
DLCS 12/13	Flexi-Cath Straight Dual Lumen Catheter 12 F, 13 cm	1
DLCS 12/16	Flexi-Cath Straight Dual Lumen Catheter 12 F, 16 cm	1
DLCS-K 12/13	Flexi-Cath Straight Dual Lumen Catheter Kit 12 F, 13 cm	1
DLCS-K 12/16	Flexi-Cath Straight Dual Lumen Catheter Kit 12 F, 16 cm	1

Domestic:

Surgimedik Healthcare India Pvt. Ltd.
Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.
Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404
domsales@surgimedik.com

International:

Surgimedik Healthcare,
223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA.
sales@surgimedik.com



SURGIMEDIK

Patient's life is our life