


# Flexi Cath T

3

## Triple Lumen Dialysis Catheter

- Ideal for extracorporeal therapies: dialysis, hemofiltration, apheresis and plasmapheresis.
- Third lumen gives greater flexibility for fluid administration, vasoactive drug therapy, nutrition and blood sampling.
- Integral soft tip reduces occurrence of stenosis and erosion.
- MRI compatible.



### Catheter Contents :

1 ea Flexi-Cath-T, 1 ea Stepper Dilator 12 F, 2 ea Injection Caps.

### Kit Contents :

1 ea Flexi-Cath-T, 1 ea Stepper Dilator 12 F, 2 ea Injection Caps.

1 ea Introducer Needle 18/7, 1 ea Guidewire JFS 80-38, 1 ea 5ml Syringe, 1 ea Scalpel.

Order No.	Description	Qty. / Box
TLCS 12/16	Flexi-Cath-T Tripple Lumen Catheter 12 F, 16 cm	1
TLCS 12/21	Flexi-Cath-T Tripple Lumen Catheter 12 F, 21 cm	1
TLCS-K 12/16	Flexi-Cath-T Tripple Lumen Catheter Kit 12 F, 16 cm	1
TLCS-K 12/21	Flexi-Cath-T Tripple Lumen Catheter Kit 12 F, 21 cm	1

### **Domestic:**

Surgimedik Healthcare India Pvt. Ltd.  
Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.  
Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404  
domsales@surgimedik.com

### **International:**

Surgimedik Healthcare,  
223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA.  
sales@surgimedik.com



**SURGIMEDIK**

Patient's life is our life