

Femo Cath

4


Single Lumen Femoral Catheter

- Material is ideal for femoral punctures.
- Available in single lumen configurations.
- Useful in over-the-wire and over-the-needle insertions.
- MRI compatible.


Luer Lock



Pinch Clamp



Soft Pliable Catheter



Catheter Contents :

1 ea Femo-Cath, 1 ea Luer Cap.

Kit Contents :

1 ea Femoral Catheter with Cap, 1 ea Guidewire, 1 ea Introducer Needle 18/7.

Order No.	Description	Qty. / Box
MFC	Femo-Cath Single Lumen Catheter 14 F, 13 cm	20
MFC-K	Femo-Cath Single Lumen Catheter Kit 14 F, 13 cm	5

Domestic:

Surgimedik Healthcare India Pvt. Ltd.
Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.
Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404
domsales@surgimedik.com

International:

Surgimedik Healthcare,
223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA.
sales@surgimedik.com



SURGIMEDIK

Patient's life is our life