


Ureteral Dilator Set

Ureteral Dilator is used for step by step dilatation of the intramural tunnel and ureter and/or for stone manipulation.



- Designed to be used over a 0.038" Guidewire
- 60 cm length with luer adaptor to permit injection

QTY. / Box : 5

Length (cm)	Description	Order no.
60	Ureteral Dilator 6F	UD 6
	Ureteral Dilator 8F	UD 8
	Ureteral Dilator 10F	UD 10
	Ureteral Dilator 12F	UD 12
	Ureteral Dilator 14F	UD 14
	Ureteral Dilator 16F	UD 16

QTY. / Box : 1

Length (cm)	Description	Order no.
60	Ureteral Dilator 6F, 8F, 10F, 12F, 14F, 16F	UDS 6/16-60

Nottingham One Step Dilator



- One Step Dilatation of the intramural portion and ureter
- May also be used to inject using female luer adaptor at proximal end
- Designed to be used over 0.038" Guidewire
- 4cm taper length, overall length 50 cm

QTY. / Box : 5

Taper From (Fr)	Taper To (Fr)	Order no.
6	10	NOT 6/10
	12	NOT 6/12

Domestic:

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.
Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

International:

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com