


PRO ECHO VIEW

Embryo Transfer Catheter



Used for ultrasound guided transfer of embryo to the uterus.

- Soft tip inner catheter
- MEA & LAL tested batchwise
- Marking on inner steel pipe & outer sheath
- Flexible transfer catheter
- Visible under ultrasound
- Soft obturator & firm stylet available separately



Soft Obturator



Stylet



Order No. (Catheter)	Order No. (Catheter with Soft Obturator)	Order No. (Catheter with Stylet)	Order No. (Soft Obturator)	Order No. (Stylet)	Qty./Box
PE-ETC	PE-ETC-SB	PE-ETC-STY	SB-PE-ETC	STY-PE-ETC	10