



# Pediatric Ureteral Stents




The Pediatric Ureteral Stent is intended to temporarily facilitate drainage of urine from kidney to bladder by trained physician.


- Extremely low encrustation tendency
- Soft polyurethane material ensures patient comfort
- Sizes to fit pediatric needs
- Contact us for other sizes



Pelvis Side  
Open Tip



Pelvis Side  
Closed Tip



Centimeter Marking  
on Stents

Stent + Pusher

Qty. / Box : 10 or 25

Open Tip Order No.	Closed Tip Order No.	Fr	Length (cm)	Compatible Guidewire
POS 3/10	PS 3/10	3	10	0.018"
POS 3/12	PS 3/12	3	12	0.018"
POS 3/14	PS 3/14	3	14	0.018"
POS 3/16	PS 3/16	3	16	0.018"
POS 4/10	PS 4/10	4	10	0.025"
POS 4/12	PS 4/12	4	12	0.025"
POS 4/14	PS 4/14	4	14	0.025"
POS 4/16	PS 4/16	4	16	0.025"

\*\* Set available with guidewire, suture, clamp etc.  
# Hydrophilic coating on stent as per requirement.

**Surgimedik  
Healthcare**

Mobile : +91-771007 9206  
Telephone : +91-22-28609401

Domestic : domsales@surgimedik.com  
International : sales@surgimedik.com

[www.surgimedik.com](http://www.surgimedik.com)