



# Ureteral Stents



The Ureteral Stent is intended to temporarily facilitate the drainage of urine from kidney to bladder by a trained physician.




- Extremely low encrustation tendency
- Soft polyurethane material ensures patient comfort
- Wide lumen for superior drainage
- Highly radiopaque for easier visualization
- 4.7 Fr Stents with long pusher for use via Ureteroscope
- Contact us for other sizes



Pelvis Side  
Open Tip



Pelvis Side  
Closed Tip



Centimeter Marking  
on Stents

Stent + Pusher

Qty. / Box : 10 or 25

Open Tip Order Number	Closed Tip Order Number	Intraoperative Tip Order Number	Fr	Length (cm)	Compatible Guidewire
SOT 4.7/16	SS 4.7/16	SSI 4.7/16	4.7	16	0.035"
SOT 4.7/24	SS 4.7/24	SSI 4.7/24	4.7	24	0.035"
SOT 4.7/26	SS 4.7/26	SSI 4.7/26	4.7	26	0.035"
SOT 4.7/28	SS 4.7/28	SSI 4.7/28	4.7	28	0.035"
SOT 4.7/30	SS 4.7/30	SSI 4.7/30	4.7	28	0.035"
SOT 5/16	SS 5/16	SSI 5/16	5	16	0.035"
SOT 5/24	SS 5/24	SSI 5/16	5	24	0.035"
SOT 5/26	SS 5/26	SSI 5/26	5	26	0.035"
SOT 5/28	SS 5/28	SSI 5/28	5	28	0.035"
SOT 5/30	SS 5/30	SSI 5/30	5	30	0.035"
SOT 6/16	SS 6/16	SSI 6/16	6	16	0.038"
SOT 6/24	SS 6/24	SSI 6/16	6	24	0.038"
SOT 6/26	SS 6/26	SSI 6/26	6	26	0.038"
SOT 6/28	SS 6/28	SSI 6/28	6	28	0.038"
SOT 6/30	SS 6/30	SSI 6/30	6	30	0.038"
SOT 7/16	SS 7/16	SSI 7/16	7	16	0.038"
SOT 7/24	SS 7/24	SSI 7/16	7	24	0.038"
SOT 7/26	SS 7/26	SSI 7/26	7	26	0.038"
SOT 7/28	SS 7/28	SSI 7/28	7	28	0.038"
SOT 7/30	SS 7/30	SSI 7/30	7	30	0.038"

\*\* Set available with guidewire, suture, clamp etc.

# Hydrophilic coating on stent as per requirement.

**Surgimedik  
Healthcare**

Mobile : +91-771007 9206  
Telephone : +91-22-28609401

Domestic : domsales@surgimedik.com  
International : sales@surgimedik.com

[www.surgimedik.com](http://www.surgimedik.com)