

# EAGLE-M


## Manual Single Lumen Ovum Pickup Needles



Used for laparoscopic or ultrasound guided transvaginal harvesting of oocyte from follicles along with flushing of follicles.

- When used with ultrasonic equipment, echo tip enhances visualization.
- Designed handle for greater control during ovum retrieval.
- Luer Lock connector for connecting a syringe
- MEA and Lal Tested
- Amount of suction applied is manually controlled.
- Available in 17

Highly echogenic tip



Order No.	Gauge	Length	Qty./Box
E-OPN-SL17/35-M	17	35 cm	10